

MADE IN EUROPE



This product has been developed and manufactured using both raw materials and with production in Europe.

Choose quality, save money, and protect the environment.

(UK) User Manual

02_20251215_UK

Sliding Material Orange



Sliding Material Orange



Article no.	HMI no.	Name - Size	UDI-DI/GTIN
MC 001-1005	811944	Sliding Material Orange - Roll 150 m	05744001970034
MC 001-1006	39037	Sliding Material Orange - Roll 10 m	05744001970041
MC 001-1009	-	Sliding Material Orange - 10 m pack	05744001970072
MC 001-1012	-	Sliding Material Orange - Roll 20 m	05744001973837
MC 001-1013	-	Sliding Material Orange - Roll 50 m	05744001973844
MC 001-1007	85291	Sliding Material Pro - Roll 10 m	05744001970058
MC 001-1008	85292	Sliding Material Pro - Roll 100 m	05744001970065

Regulation MDR 2017/745 for medical devices class 1



ISO 12 31 03 21

The Original Nordic Human Handling Concept



Developed & Manufactured by:



Thank you for choosing a product from master care A/S.

We recommend reading the user manual thoroughly before using the product (-s).

The information in this user manual applies to the shown products provided that the products have not been modified or customized and that the products have been used as recommended in this user manual.

If the customer changes, repairs or uses the product in other combinations than recommended by master care A/S the CE-approval no longer applies and the guarantee becomes void.

In case of doubt, contact master care A/S or our distributor in your area.

MDR (EU) 2017/745:

Products in this user manual are CE-approved and follows the harmonized standards. This means that the products are in compliance with the essential health and safety directives. The product(-s) comply with the standards in the Medical Device Regulation MDR (EU) 2017/745 – Class 1.

The intended use of the product/-s:

The product is to be used bilayered to reduce friction for users with decreased function during e.g. transfer, dressing or training. The product can be used to reduce the risk of shear injuries.

Correct use of the product:

Place the Nylon underneath the points/support areas where friction is to be reduced. Always use the product double layered to achieve maximum functionality and reduction of friction. Always remove the product under the user after use to avoid fall injuries and risk of skin overheating.

Warning:

Always do a risk assessment before using the product especially if used e.g., when sitting on the edge of the bed since this poses a risk of the user sliding/falling on to the floor. Always make sure the user has sufficient support underneath their feet when sitting if a nylon

piece is placed underneath their bottom. Never leave the product on the floor since this poses a risk of falling.

Risk assessment:

Always make a risk assessment and a functional assessment of the user prior to using the product. This way using the product for transferring or moving a user, is ergonomic and legally sound.

Knowledge of the product/transferring:

It is important that the user of the product has been properly instructed in the use of the product and has read the user manual thoroughly prior to using the product.

master care A/S and our distributors provide education in the use of products and on our website all user manuals can be found.

Life span:

The expected life span of products with nylon sliding material with daily wear and tear and washing averagely is 1-3 years. The products life span depends on the use, users load, washing and drying temperatures and how often it is washed and dried.

Make visual inspections frequently and preferably before every use. Products with insufficient friction reduction, defect seams, broken material or other deficiencies must be discarded. This way you will always keep a sound and healthy work environment for both helper and user.

Guarantee:

master care A/S provide a 2-year guarantee covering material defects and manufacturing defects if the product has been used correctly.

Unintended events and injuries when misusing the products:

Unintended events must always be reported to the appropriate National Authority and to master care A/S or the distributor in the current country.

master care A/S and our distributors take no responsibility for personal injuries or product damages if the product(-s) have been used in contravention of the by the manufacturer recommended intended purpose or proper use of the product.



User Manual - Sliding Material Orange and Pro

Table of Contents

Turning and repositioning in bed:	page 2
Getting further up in bed, more active user:	page 2
Getting further up in bed, less active user:	page 2
Move the user higher up ontop of the Master Turner:.....	page 3
Transferring from chair and bed:	page 3
Repositioning with Multi Sheet:.....	page 3
Technical Specifikations	page 4




master

User Manual - Sliding Material Orange and Pro


Turning and repositioning in bed:



1. Cut a sufficient piece of Nylon from the roll. Fold the Nylon piece and place the fold under the user, along the spine or under the draw sheet towards the side of the transfer.



2. Bend the user's legs. Place the user's hands on the opposite shoulder or on the pillow. Return the user to supine lying on the Nylon piece and slide/reposition the user to side lying.



3. Remove the Nylon piece by pulling diagonally from the lower part. If you feel a resistance just pull from the opposite lower corner.

Getting further up in bed, more active user:



4. Fold the Nylon piece and place the hand in the fold. From above slide the Nylon under the user's shoulders – first on one side and afterwards on the other side.



5. Ask the user to bend the knees and place feet on the sheet. By pushing the feet down and lifting the pelvis from the sheet the user slides the upper body upwards. If necessary, support the user's feet.



6. Remove the Nylon piece by pulling on the corner of the lower part of the material.

Getting further up in bed, less active user:



7. Turn the user to side lying position. Fold the Nylon piece and place it with the fold under the pelvis. It must cover the pressure points under the head, shoulders and pelvis as the user is turned back to supine lying. If necessary, straighten out the Nylon before the transfer.




8. Ask the user to bend the knees and place the feet on the sheet. Elevate the foot of the bed to use the force of gravity. If necessary, support the user's feet.
Alternatively: Place the Nylon piece under the sheet and use this as a draw sheet during the transfer.



9. Remove the Nylon piece by pulling on the corner of the lower part of the material. If you feel a resistance just pull from the opposite lower corner.

User Manual - Sliding Material Orange and Pro

Move the user higher up onto of the Master Turner:



10. Turn the user to side lying position. Fold the Nylon piece and place it with the fold under the pelvis. It must cover the pressure points under the head, shoulders and pelvis as the user is turned back to supine lying.



11. Fix the Master Turner towards the Mattress Cover. Ask user to bend the knees and place the feet on the Mattress Cover for the user to move independently. Alternatively stabilize the user's feet onto the Master Turner/keep the legs straight, elevate the foot end of the bed and push shin/feet to transfer upwards.



12. Remove the Nylon piece by pulling on the corner of the lower part of the material. If you feel a resistance just pull from the opposite lower corner.

Transferring from chair and bed:



13. Ask the user to lean away from the bed. Fold the Nylon piece and slide it under the elevated part of the pelvis.



14. Place the Master Glide sliding board under the folded Nylon piece to build a bridge and make the transferring easier. Remember to make use of the user's resources during the transfer.



15. Remove the sliding board and afterwards remove the Nylon piece by pulling on the corner of the lower part of the material. If you feel a resistance just pull from the opposite lower corner.

Repositioning with Multi Sheet:



16. Turn the user to side lying position using the Multi Sheet - if necessary, connected to the ceiling hoist. Fold the Nylon piece and place it with the fold towards the spine. Make sure it covers all the pressure points i.e., pelvis, shoulders, head.



17. Connect the Multi Sheet to the hoist on the opposite side. Elevate and slide the user into position using the Multi Sheet as a draw sheet.



18. Remove the Nylon piece by pulling on the corner of the lower part of the material. If you feel a resistance just pull from the opposite lower corner.

User Manual - Sliding Material Orange and Pro

Technical Specifications

Materials:

Single side coated Nylon

ISO certifications:

DS/EN ISO 12952-1:2015
Textiles – Assessing the ignitability of all bedding items - Part 1: Ignition source: smouldering cigarette.

Article no. - Product Name:

Size:

Article no. - Product Name:	Size:
MC 001-1005 Sliding Material Orange - Roll 150 m	L 150 x W 1,6 m
MC 001-1006 Sliding Material Orange - Roll 10 m	L 10 x W 1,6 m
MC 001-1009 Sliding Material Orange - 10 m pack	L 10 x W 1,6 m
MC 001-1012 Sliding Material Orange - Roll 20 m	L 20 x W 1,6 m
MC 001-1013 Sliding Material Orange - Roll 50 m	L 50 x W 1,6 m
MC 001-1007 Sliding Material Pro - Roll 10 m	L 10 x W 1,6 m
MC 001-1008 Sliding Material Pro - Roll 100 m	L 100 x W 1,6 m

WARNINGS



Never leave the user on the product when lying on the edge of the bed.



Never leave the product on the floor.



Make sure that all relevant brakes are activated before transferring the user.



CE-certification:

This product is CE-approved and follows the harmonized standards. This means that the products are in compliance with the essential health and safety directives. The product(-s) comply with the standards in the Medical Device Regulation (EU) 2017/745 – class 1.



User Weight:

No maximum weight limit.



Ironing:

Do not iron.



Washing:

This product can be washed in a washing machine at max. 85 °C. Gentle washing and dry tumbling significantly prolongs the durability.



Dry cleaning:

Do not dry clean.



Detergents pH level:

This product tolerates detergents / disinfectants with a pH level between 5-9.



Bleach:

Do not bleach.



Fabric softener:

Never use fabric softener as this will reduce the sliding properties significantly.



Maintenance:

Make visual inspections frequently and preferably before every use. Pay attention to torn fabric, seams and straps. Products with high friction or other deficiencies must be discarded.



Dry Tumbling:

This product may only be dry tumbled gently and at max. 40 °C. Airdrying prolongs the durability significantly.



User Manual:

Always read the User Manual thoroughly before using the product.



Surface cleaning:

The surface of this product can be cleaned with a wet cloth and soap.



Disposal:

Always make sure to follow local guidelines when disposing or discarding the product.



Disinfectants:

This product tolerates alcoholic disinfectants with a pH level between 5-9.

For more information, videos Quick Guide etc., visit www.master-care.dk

Find us on:

