

# MADE IN EUROPE



This product has been developed and manufactured using both raw materials and with production in Europe.

Choose quality, save money, and protect the environment.

# (UK) User Manual

02\_20221209\_UK

## Nylon Master



MC 001-1101



MC 001-1106



MC 001-1204




MC 001-1206



MC 001-1220



MC 001-1271



Article no.	HMI no.	Name - Size	UDI-DI/GTIN
<input type="radio"/> MC 001-1101 <sup>1</sup>	34304	Nylon Master Tube - XS	5744001970171
<input type="radio"/> MC 001-1102 <sup>1</sup>	34305	Nylon Master Tube - S	5744001970188
<input type="radio"/> MC 001-1103 <sup>1</sup>	34306	Nylon Master Tube - M	5744001970195
<input type="radio"/> MC 001-1104 <sup>2</sup>	34307	Nylon Master Tube - L	5744001970201
<input type="radio"/> MC 001-1106 <sup>1</sup>	81812	Nylon Master Tube - Basic	5744001973097
<input type="radio"/> MC 001-1201 <sup>1</sup>	34300	Nylon Master Open - XS	5744001973103
<input type="radio"/> MC 001-1202 <sup>1</sup>	34301	Nylon Master Open - S	5744001970232
<input type="radio"/> MC 001-1203 <sup>1</sup>	34302	Nylon Master Open - M	5744001970249
<input type="radio"/> MC 001-1204 <sup>2</sup>	34303	Nylon Master Open - L	5744001970256
<input type="radio"/> MC 001-1206 <sup>2</sup>	101812	Nylon Master Open - Overlock	5744001970263
<input type="radio"/> MC 001-1210 <sup>1</sup>	33232	Slide Master Nylon Padded - S	5744001973127
<input type="radio"/> MC 001-1220 <sup>1</sup>	33233	Slide Master Nylon Padded - M	5744001970270
<input type="radio"/> MC 001-1271 <sup>2</sup>	39025	Slide Master Nylon Padded - Wide	5744001973134

Regulation MDR 2017/745 for medical devices class 1



<sup>1</sup> ISO 12 31 03 20

<sup>2</sup> ISO 12 31 03 21

The Original Nordic  
Human Handling Concept



Developed &  
Manufactured by:



# Tak for at du har valgt et produkt fra master care A/S.

We recommend reading the user manual thoroughly before using the product (-s).

The information in this user manual applies to the shown products provided that the products have not been modified or customized and that the products have been used as recommended in this user manual.

If the costumer changes, repairs or uses the product in other combinations than recommended by master care A/S the CE-approval no longer applies and the guarantee becomes void.

In case of doubt, contact master care A/S or our distributor in your area.

## **MDR (EU) 2017/745:**

Products in this user manual are CE-approved and follows the harmonized standards. This means that the products are in compliance with the essential health and safety directives. The product(-s) comply with the standards in the Medical Device Regulation MDR (EU) 2017/745 – Class 1.

## **The intended use of the product/-s:**

The product is to be used bilayed to reduce friction for users with decreased function during e.g. transfer, dressing or training. The product can be used to reduce the risk of shear injuries.

## **Correct use of the product:**

Place the Nylon underneath the points/support areas where friction is to be reduced. Always use the product double layered to achieve maximum functionality and reduction of friction. Always remove the product under the user after use to avoid fall injuries and risk of skin overheating.

## **Warning:**

Always do a risk assessment before using the product especially if used e.g., when sitting on the edge of the bed since this poses a risk of the user sliding/falling on to the floor. Always make sure the user has sufficient support underneath their feet when sitting if a nylon piece is placed underneath their bottom.

Never leave the product on the floor since this poses a risk of falling.

## **Risk assessment:**

Always make a risk assessment and a functional assessment of the user prior to using the product. This way using the product for transferring or moving a user, is ergonomic and legally sound.

## **Knowledge of the product/transferring:**

It is important that the user of the product has been properly instructed in the use of the product and has read the user manual thoroughly prior to using the product.

master care A/S and our distributors provide education in the use of products and on our website all user manuals can be found.

## **Life span:**

The expected life span of products with nylon sliding material with daily wear and tear and washing averagely is 1-3 years. The products life span depends on the use, users load, washing and drying temperatures and how often it is washed and dried.

Make visual inspections frequently and preferably before every use. Products with insufficient friction reduction, defect seams, broken material or other deficiencies must be discarded. This way you will always keep a sound and healthy work environment for both helper and user.

## **Guarantee:**

master care A/S provide a 2-year guarantee covering material defects and manufacturing defects if the product has been used correctly.

## **Unintended events and injuries when misusing the products:**

Unintended events must always be reported to the appropriate National Authority and to master care A/S or the distributor in the current country.

master care A/S and our distributors take no responsibility for personal injuries or product damages if the product(-s) have been used in contravention of the by the manufacturer recommended intended purpose or proper use of the product.



# User Manual - Nylon Master

## Table of Contents

Slide Master – Positioning the Nylon Master	.....	page 4
Turning in the bed with the Nylon Master/Slide Master	.....	page 4
Slide Master - Moving higher up in the bed with the Nylon Master	.....	page 4
Using Nylon Master to bend the legs	.....	page 5
Move the pillow using Nylon Master	.....	page 5
Remove the Nylon Master/Slide Master	.....	page 5
Technical Specifications	.....	page 6



# User Manual - Nylon Master

## Slide Master – Positioning the Nylon Master:



1. The patient bends his/her leg and grabs the bed rail. The Nylon Master/Slide Master is folded and inserted under the patient's pressure points.



2. The patient grips the opposite bed rail and the helper takes hold of the Nylon Master/Slide Master and straightens it out.



3. The patient returns to a supine position and the Nylon Master/Slide Master is now positioned correctly under the patient. Friction under the patient's pressure points has been reduced.

## Turning in the bed with the Nylon Master/Slide Master:



4. The Nylon Master/Slide Master is placed under the patient's pressure points.



5. Turn the user to the side and pull Nylon Master/Sling Master lightly towards you.



6. The user bends his legs. The helper places his hands on the user's hip and knee and turns the user to side position.

## Slide Master - Moving higher up in the bed with the Nylon Master:



7. The Nylon Master/Slide Master is placed under the patient's shoulder and pelvis.



8. The patient bends his/her legs and, pushing off, lifts his/her pelvis. The foot end may be elevated.



# User Manual - Nylon Master

## Using Nylon Master to bend the legs:



9. Place the Nylon Master bilayed under the user's lower legs and feet, with the fold towards the body. Use the knee bend or the free space over the heel to place the Nylon Master. Support the foot and knee when bending the leg while sliding on top of the Nylon Master.

## Move the pillow using Nylon Master:



10. Ideal for user's who pushes the head into the pillow or cannot move the head independently. Place the Nylon Master between the user's head and the pillow. Support the user's head while moving the pillow.

## Remove the Nylon Master/Slide Master:



11. Push the Nylon Master / Slide Master under the lower back and bottom on one side of the user.



12. On the other side the helper takes hold of the bottom layer of the Nylon Master/Slide Master.



13. The helper pulls the Nylon Master/Slide Master out using weight transfer, such that the double slide is utilized.

# User Manual - Nylon Master

## Technical Specifications

### Materials:

**Nylon Master:** Nylon

**Slide Master:** Nylon, polyester

Article no. - Product Name	Size:	Max. User Weight:
MC 001-1101 Nylon Master Tube - XS	L 30 x W 60 cm	180 kg
MC 001-1102 Nylon Master Tube - S	L 60 x W 60 cm	180 kg
MC 001-1103 Nylon Master Tube - M	L 80 x W 80 cm	180 kg
MC 001-1104 Nylon Master Tube - L	L 100 x W 80 cm	180 kg
MC 001-1106 Nylon Master Tube - Basic	L 75 x W 70 cm	180 kg
MC 001-1201 Nylon Master Open - XS	L 30 x W 120 cm	180 kg
MC 001-1202 Nylon Master Open - S	L 60 x W 120 cm	180 kg
MC 001-1203 Nylon Master Open - M	L 80 x W 160 cm	180 kg
MC 001-1204 Nylon Master Open - L	L 100 x W 160 cm	180 kg
MC 001-1206 Nylon Master Open - Overlock	L 80 x W 160 cm	180 kg
MC 001-1210 Slide Master Nylon Padded - S	L 30 x W 60 cm	180 kg
MC 001-1220 Slide Master Nylon Padded - M	L 50 x W 60 cm	180 kg
MC 001-1271 Slide Master Nylon Padded - Wide	L 60 x W 100 cm	180 kg

## WARNINGS



Never leave the user on the product when lying on the edge of the bed.



Never leave the product on the floor.



Make sure that all relevant brakes are activated before transferring the user.



### CE-certification:

This product is CE-approved and follows the harmonized standards. This means that the products are in compliance with the essential health and safety directives. The product(-s) comply with the standards in the Medical Device Regulation (EU) 2017/745 – class 1.



### User Weight:

Recommended max. user weight: 180 kg



### Ironing:

Do not iron.



### Washing:

This product can be washed in a washing machine at max 85 °C. Gentle washing and dry tumbling significantly prolongs the durability.



### Dry cleaning:

Do not dry clean.



### Detergents pH level:

This product tolerates detergents / disinfectants with a pH level between 5-9.



### Bleach:

Do not bleach.



### Fabric softener:

Never use fabric softener as this will reduce the sliding properties significantly.



### Maintenance:

Make visual inspections frequently and preferably before every use. Pay attention to torn fabric and seams etc. Products with high friction or other deficiencies must be discarded.



### Dry Tumbling:

This product may only be dry tumbled gently and at max. 40 °C. Airdrying prolongs the durability significantly.



### User Manual:

Always read the User Manual thoroughly before using the product.



### Surface cleaning:

The surface of this product can be cleaned with a wet cloth and soap.



### Disposal:

Always make sure to follow local guidelines when disposing or discarding the product.



### Disinfectants:

This product tolerates alcoholic disinfectants with a pH level between 5-9.

Find us on:

