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○ **MC 001-1951** [ HMI: 81134 ]  
Mattress Cover Extender -  
**Large**

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

**UK** Place cover at the head of the mattress with the label "This side up and towards head end." Hold the cover firmly against the mattress and fold the mattress up and down against the foot of the bed.

Close the ends of the mattress cover with the attached Velcro-like strips and tighten. Fold the mattress from the foot to the head and close the ends in the same way


The mattress cover is now fitted. Tighten and check at least once a day.

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
**UK** If you are fitting the mattress cover on a pressure relieving mattress, you can purchase 4 mattress cover extensions

If there is a mattress topper on top of a regular mattress, fit the mattress cover around both mattresses

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

DK 2 folder derefter den øverste strop henover den nederste strop

UK 1 holds on to the Master Turner, while 2 loosens the two Velcro-like strips on the sliding flap at the head

UK 2 folds the bottom strap and pulls it up towards the user's neck and thereby away under the user's back

UK 2 then folds the top strap across the bottom strap

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

UK 2 places his/her palm on top of the strap and tucks it in under the user's neck and over towards the opposite side

UK 1 takes the strap and now pulls the sliding flap out with both hands

UK The Master Turner shuts it against the anti-slip lining. The user's body is now secured against the surface

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
**UK** 2 takes hold of the Master Turner and tightens it around the user's body


1 places the sliding flap up against the Master turner and releases it. 1 places his/her palms on the sliding flap and tucks it in between the mattress cover and the Master Turner

1 first tucks the sliding flap under the user in the middle and smooths it out towards the sides. The straps are tucked on top of the sliding flap and under the user's natural cavities (neck, lower back, back of the knee and over the heel)

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**UK** Cross the user's legs. Guide hand to opposite shoulder. Push the pillow and shoulder forward on the side, to which the user is being turned

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UK

- Easy fittin
- The slide flaps can fix the turn sheet as need
- Works with the bed mattress
- Does not diminish the effects of alternating pressure mattresses
- 10 cm non-slip coating along the edge on both sides
- Unique sliding properties - HLF (Hyper Low Friction)
- Prevents pressure ulcers
- Extremely durable

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**MC 001-1732** [ HMI: 49368 ]  
Master Turner Soft - **Large**

Størrelse / Dimensions / Größe:  
L 215 x B/W 250

**MC 001-2540** [ HMI: 49370 ]  
Master Turner Two In One Soft - **Large**

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
Størrelse / Dimensions / Größe:  
L 47 x B/W 4 cm

**MC 001-2500** [ HMI: 49873 ]  
Master Turner DUO Two In One Soft Bund / Bottom  
/ Unterteil - **Large**

Størrelse / Dimensions / Größe:  
L 160 x B/W 145 cm

**MC 001-2501** [ HMI: 49872 ]  
Master Turner DUO Two In One Soft Top / Top /  
Oberteil - **Large**

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**MC 001-1951** [ HMI: 81134 ]  
Mattress cover extender - **Large**

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Washing temperature: **40 - 85°C**

Tumble dry: **max 60°C**

Airdrying: extends durability significantl

Do not use softener.

Use detergent without bleach and chlorine-containing compounds.

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### WARNING

There is a risk of the user slipping over the edge of the bed

Never leave the user lying on master **care's** products with lowered safety rails

Never leave the user with master **care's** products attached to the lift or with lowered safety rails

### IMPORTANT

- These products must not be used if they show critical wear or are otherwise damaged
- Always use the products in accordance with this manual
- Always check the stitching on these products - both before and after use
- If in doubt about the proper use of the products, please contact master **care**/your dealer immediately

- master **care** assumes no liability for accidents resulting in injuries if these products are not used in accordance with this manual
- Trading with master **care** is conducted in accordance with the company's at all time applicable Trade and Business Conditions
- **Velcro straps must be closed by washing and drying**
- **During wash and drying, the washing straps MUST be applied to the velcro on the toppiece**

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These products comply with the standards for Class 1 products in the Medical Device Regulation (EU) 2017/745.

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